

Product Overview



- Adjustable height from 71 cm to 111 cm (27.95 inches to 43.70 inches)
- 2 liftable columns with a weight capacity of 160 Kg in total
- Worksurface from W 100 to 140 x D80 cm (W 39.4 to 55, Depth 31.5 inches)
- Feedtroughs openings (3holes) for cables and tubing
- Worksurface can easily be lifted up or down with retractable control panel LCD (4 positions storable)
- Recommended use in the Liquid Chromatography field
- Caster wheels allow to move the IonLC desk easily
- Four lockable castor wheels
- AC 110 60 Hz or 220 V 50 Hz



Retractable Control Module



DR2 Dedicated Drawers




IECKIT - IEC Power Strip + cables



ARM1 Information Package



PC1 UC Holder



Upper position



Lower position



Description	Part-number
IonLC-Desk - 100 cm width x80 cm depth IonLC-Desk - 39.37 inch width x 31.50 inch depth (220V / 50Hz)	IonLC-Desk-10080-EU
IonLC-Desk - 120 cm width x80 cm depth IonLC-Desk - 47.24 inch width x 31.50 inch depth (220V / 50Hz)	IonLC-Desk-12080_EU
IonLC-Desk - 140 cm width x 80 cm depth IonLC-Desk - 55.12 inch width x 31.50 inch depth (220V / 50Hz)	IonLC-Desk-14080_EU
IonLC-Desk - 100 cm width x80 cm IonLC-Desk - 39.37 inch width x 31.50 inch depth (110V / 60Hz)	IonLC-Desk-10080_US
IonLC-Desk - 120 cm width x80 cm depth IonLC-Desk - 47.24 inch width x 31.50 inch depth (110V / 60Hz)	IonLC-Desk-12080_US
IonLC-Desk - 140 cm width x 80 cm depth IonLC-Desk - 55.12 inch width x 31.50 inch depth (110V / 60Hz)	IonLC-Desk-14080_US
Optional Accessories	Part-number
CPU Holder	PC1
IEC power strip (6 outlets, and the x 6 cables are provided)	IECKIT
Solvent waste container 30 cm x 20 cm x 20 cm	SOLV
Set of 2 pull out drawers	DR2
Optional Information Control Package (Flat screen monitor arm, keyboard & mouse platform)	ARM1